


PRODUCT: PBS (10X), pH 7.4
Phosphate Buffered Saline

CATALOG#: 119-069-491 **Size:** 4L

LOT#: 723726

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Result</u>
pH:	For information only	7.415 (diluted to 1x)
OSMOLALITY:	For information only	296 (diluted to 1x) Osm/kg of water
ENDOTOXIN:	Test & Report	<0.01 EU/ml (diluted to 1x)
MYCOPLASMA:	Negative	pass (diluted to 1x)
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	3/2020	
EXPIRATION DATE:	3/2022	



Quality Assurance

4/15/2020

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.