

**PRODUCT:** DMEM  
Dulbecco's Modified Eagle's Medium  
With L-glutamine

**CATALOG#:** 112-014-101      **Size:** 500ml


**LOT#:** 723806

---

**CERTIFICATE OF ANALYSIS**

---

	<u>Range</u>	<u>Result</u>
<b>pH:</b>	6.8 – 7.4	7.152
<b>OSMOLALITY:</b>	318 - 352 mOsm/kg of water	339
<b>ENDOTOXIN:</b>	Test & Report	<0.020 EU/ml
<b>MYCOPLASMA:</b>	Negative	pass
<b>APPEARANCE:</b>	Pink to a reddish-orange colored fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
<b>CELL GROWTH:</b>	Grows mostly monolayer cells. Cells have been propagated for two consecutive passage with an average 3.0/4.0 fold increase in four days.	pass
<b>RECOMMENDED STORAGE CONDITIONS:</b>	2 – 8°C	
<b>MANUFACTURED DATE:</b>	6/2020	
<b>EXPIRATION DATE:</b>	6/2021	

  
Quality Assurance

6/29/2020

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**