


PRODUCT: NORMAL SALINE
With 9.0 g/Liter Sodium Chloride

CATALOG#: 114-055-101 Size: 500ml

LOT#: 723449

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
OSMOLALITY:	276 - 305	286
APPEARANCE:	Clear colorless fluid.	Pass
ENDOTOXIN:	Test & Report	<0.028 EU/ml
MYCOPLASMA:	Negative	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	9/2019	
EXPIRATION DATE:	9/2021	


Quality Assurance

10/2/2019

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**