

**PRODUCT:** TBS (10X), pH 7.4  
Tris Buffered Saline

**CATALOG#:** 351-086-151      Size: 10L

**LOT#:** 724000

---

**CERTIFICATE OF ANALYSIS**

---

|                                            | <u>Specification</u>    | <u>Results</u> |
|--------------------------------------------|-------------------------|----------------|
| <b>pH:</b>                                 | 7.3 – 7.5               | 7.362          |
| <b>CONDUCTIVITY:</b>                       | For information only    | 116.0 mS/cm    |
| <b>APPEARANCE:</b>                         | Clear, colorless liquid | pass           |
| <b>DNASE:</b>                              | Not detectable          | pass           |
| <b>RNASE:</b>                              | Not detectable          | pass           |
| <b>PROTEASE:</b>                           | Not detectable          | pass           |
| <b>RECOMMENDED<br/>STORAGE CONDITIONS:</b> | 15 – 30°C               |                |
| <b>MANUFACTURED DATE:</b>                  | 10/2020                 |                |
| <b>EXPIRATION DATE:</b>                    | 10/2023                 |                |

*Ta'Shemma Jackson*

10/28/2020

---

Quality Assurance

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**