

PRODUCT: RPMI-1640 MEDIUM with L-Glutamine

CATALOG#: 112-025-101 Size: 500ml

LOT#: 723976

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.2 – 7.4	7.233
OSMOLALITY:	266 - 294 mOsm/kg of water	283
ENDOTOXIN:	Test & Report	0.049 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
CELL GROWTH:	Grows mostly suspension cells. Cells have been propagated for two consecutive passages with an average 3.0-4.0 fold increase in four days.	pass

RECOMMENDED

STORAGE CONDITIONS: 2 – 8°C

MANUFACTURED DATE: 10/2020

EXPIRATION DATE: 10/2022

Ta'Shemma Jackson 10/23/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**