

PRODUCT: WATER, Cell Culture
CATALOG#: 118-162-161 Size: 20L
LOT#: 723973

CERTIFICATE OF ANALYSIS

This product is de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Result</u>
RESISTIVITY:	15 – 20 Meg Ohms/cm	17.95 Meg Ohms/cm
ENDOTOXIN:	Less than 0.03 Eu/ml	pass
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
CELL GROWTH:	Grows mostly monolayer cells. Cells have been propagated for two consecutive passages with an average 3.0/4.0 fold increase in four days.	pass

RECOMMENDED

STORAGE CONDITIONS: 15 – 30°C
MANUFACTURED DATE: 9/2020
EXPIRATION DATE: 9/2022

Michael J. Fisher

10/21/2020

Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.