

PRODUCT: DPBS without Ca & Mg
Dulbecco's Phosphate Buffered Saline

CATALOG#: 114-057-491 Size: 4L

LOT#: 723970

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.2 – 7.4	7.283
OSMOLALITY:	266 - 294 mOsm/kg of water	285
ENDOTOXIN:	Test & Report	<0.028 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	9/2020	
EXPIRATION DATE:	9/2022	

Michael J. Fisher

Quality Assurance

10/21/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**