

PRODUCT: Hypotonic Buffer (1mM TRIS buffer, pH 9.0)

CATALOG#: 351-374-136 Size: 20L

LOT#: 723967

CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
pH:	8.9 – 9.1	8.960 @ 21.6°C
CONDUCTIVITY:	Report	36 µS/cm @ 22.3°C
APPEARANCE:	Clear, colorless liquid	pass
DENSITY	Report	0.9996 grams/ml
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	9/28/2020	
EXPIRATION DATE:	9/2022	

Ta'Shemma Jackson

10/22/2020

Quality Assurance

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY