

PRODUCT: RPMI-1640 Medium with Glutamine XL
CATALOG#: 112-301-101 **Size:** 500ml
LOT#: 723956

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.0 – 7.4	7.278
OSMOLALITY:	266 - 294 mOsm/kg of water	284
ENDOTOXIN:	Test & Report	0.068 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass

**RECOMMENDED
STORAGE CONDITIONS:** 2 – 8°C
MANUFACTURED DATE: 9/2020
EXPIRATION DATE: 9/2022

Michael J. Fisher

10/09/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**