

PRODUCT: Tris-HCl (1M), pH 8.0

CATALOG #: 351-007-101

Size: 500ml

LOT #: 723947

CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
pH:	7.9 – 8.1	7.970
Conductivity:	For information only	36,100 μ S
Appearance:	Clear colorless liquid	pass
DNase:	Not detectable	pass
RNase:	Not detectable	pass
Protease:	Not detectable at 0.5M concentration pass	

RECOMMENDED

STORAGE CONDITIONS: 15 – 30°C

MANUFACTURED DATE: 9/2020

EXPIRATION DATE: 9/2023

Michael J. Fisher

9/18/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**