

PRODUCT: PBS (1X), pH 7.4
Phosphate Buffered Saline

CATALOG#: 114-058-151 Size: 10L

LOT#: 723946

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.3 – 7.5	7.378
OSMOLALITY:	285 - 315 mOsm/Kg H ₂ O	296
ENDOTOXIN:	Test & Report	0.042 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	9/2020	
EXPIRATION DATE:	9/2022	

Michael J. Fisher

10/01/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**