

PRODUCT: LB BROTH (LURIA-BERTANI, Miller)

CATALOG #: 340-004-101 **Size:** 500ml

LOT #: 723905

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
pH:	7.0 - 7.4	7.180
OSMOLALITY:	389 - 431 mOsm/Kg of water	411
APPEARANCE:	Light amber, clear to slightly opalescent liquid	pass
STERILITY:	Sterilized using saturated steam under pressure	pass
PERFORMANCE:	Passed E. coli growth performance test.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	8/2020	
EXPIRATION DATE:	8/2022	

Chanika Johnson

8/19/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**