

PRODUCT: PBS (1X), pH 7.4
Phosphate Buffered Saline

CATALOG#: 114-058-101 Size: 500ml

LOT#: 723888

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.3 – 7.5	7.381
OSMOLALITY:	285 - 315 mOsm/Kg H ₂ O	298
MYCOPLASMA:	Negative	pass
ENDOTOXIN:	Test & Report	0.031 EU/ml
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Michael J. Fisher

8/21/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**