

**PRODUCT:** WATER, Endotoxin Free (<0.005 EU/ml)  
**CATALOG#:** 118-325-101      Size: 500ml  
**LOT#:** 723886

---

### **CERTIFICATE OF ANALYSIS**

---

This product is de-ionized by reverse osmosis and sterile filtered.

	<b><u>Specifications</u></b>	<b><u>Result</u></b>
<b>RESISTIVITY:</b>	15 – 20 Meg Ohms/cm	17.95 Meg Ohms/cm
<b>ENDOTOXIN:</b>	Less than 0.005 Eu/ml	pass
<b>MYCOPLASMA:</b>	Negative	pass
<b>APPEARANCE:</b>	Clear colorless fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
<b>DNASE:</b>	Not detectable	pass
<b>RNASE:</b>	Not detectable	pass
<b>PROTEASE:</b>	Not detectable	pass
<b>RECOMMENDED STORAGE CONDITIONS:</b>	15 – 30°C	
<b>MANUFACTURED DATE:</b>	7/2020	
<b>EXPIRATION DATE:</b>	7/2022	

*Chanika Johnson*

8/19/2020

---

Quality Assurance

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**