


PRODUCT: DPBS (1X), pH 7.3-7.5 without Ca & Mg
CATALOG#: 114-309-151 **Size:** 10L
LOT#: 723002

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.4±0.1	7.370
OSMOLALITY:	280 ± 5% mOsm/kg of water	282
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	12/2018	
EXPIRATION DATE:	12/2020	



Quality Assurance

12/27/2018

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**