

PRODUCT: TBE Buffer 10X
Tris Borate EDTA Buffer

CATALOG#: 351-001-491 Size: 4L

LOT#: 723123


CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
pH:	8.1 – 8.5	8.359
CONDUCTIVITY:	4.3 – 4.7 mS/cm	4.71 mS/cm
APPEARANCE:	Clear, colorless liquid	pass
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass

**RECOMMENDED
STORAGE CONDITIONS:** 15 – 30°C

MANUFACTURED DATE: 4/2019

EXPIRATION DATE: 4/2022


Quality Assurance

4/15/2019

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**