


PRODUCT: ACK LYSING BUFFER
CATALOG#: 118-156-101 Size: 500ml
LOT#: 723165

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.2 – 7.6	7.520
OSMOLALITY :	276 - 305 mOsm/Kg H ₂ O	303 mOsm/Kg H ₂ O
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass

RECOMMENDED

STORAGE CONDITIONS: 15 – 30°C
MANUFACTURED DATE: 3/2019
EXPIRATION DATE: 3/2021



Quality Assurance

4/15/2019

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**