


PRODUCT: DPBS 10X
Dulbecco's Phosphate Buffered Saline
without Calcium and Magnesium

CATALOG#: 119-068-151 **Size:** 10L

LOT#: 723177

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.2 – 7.4 (diluted to 1x)	7.326
OSMOLALITY:	266 - 294 (diluted to 1x) Osm/kg of water	280
ENDOTOXIN:	Test & Report	<0.020 EU/ml
MYCOPLASMA:	Negative	pass (diluted to 1x)
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	4/2019	
EXPIRATION DATE:	4/2021	


Quality Assurance

5/7/2019

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**