


PRODUCT: DPBS (1X), pH 7.3-7.5 without Ca & Mg
CATALOG#: 114-309-151 **Size:** 10L
LOT#: 723360

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.4±0.1	7.390
OSMOLALITY:	280 ± 5% mOsm/kg of water	280
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2019	
EXPIRATION DATE:	7/2021	



Quality Assurance

8/13/2019

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**