

PRODUCT: WATER, Cell Culture
CATALOG#: 118-162-151 Size: 10L
LOT#: 724131

CERTIFICATE OF ANALYSIS

This product is de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Result</u>
RESISTIVITY:	15 – 20 Meg Ohms/cm	17.95 Meg Ohms/cm
ENDOTOXIN:	Less than 0.03 Eu/ml	pass
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	2/2021	
EXPIRATION DATE:	2/2023	

Michael J. Fisher

Quality Assurance

02/23/2021

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.