

PRODUCT: LB BROTH (LENNOX)
CATALOG#: 340-016-131 Size: 1000ml
LOT#: 723856

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
pH:	7.0 – 7.4	7.176
OSMOLALITY:	235 - 260 mOsm/Kg of water	252
APPEARANCE:	Light amber, clear to slightly opalescent liquid	pass
STERILITY:	Sterilized using saturated steam under pressure	pass
PERFORMANCE:	Passed E. coli growth performance test.	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Michael J. Fisher

7/23/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**