

PRODUCT: RPMI-1640 Medium with Glutamine XL

CATALOG#: 112-301-101 Size: 500ml

LOT#: 723855

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.0 – 7.4	7.229
OSMOLALITY:	266 - 294 mOsm/kg of water	287
ENDOTOXIN:	Test & Report	0.062 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass

RECOMMENDED

STORAGE CONDITIONS: 2 – 8°C

MANUFACTURED DATE: 7/2020

EXPIRATION DATE: 7/2022

Michael J. Fisher

8/05/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**