

PRODUCT: SUPER BROTH
CATALOG #: 340-018-131 **Size:** 1000ml
LOT #: 723844

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
pH:	6.8 – 7.2	7.092
OSMOLALITY:	For information only	529 mOsm/kg of water
APPEARANCE:	Clear, light to medium amber liquid	pass
STERILITY:	Sterilized using saturated steam under pressure	pass
PERFORMANCE:	Passed E. coli growth performance test.	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Michael J. Fisher

7/30/2020

Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.