

PRODUCT: AMPHOTERICIN B 250 µg/ml
CATALOG#: 120-097-711 Size: 4x50ml, 50ml
LOT#: 723839

CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
pH:	For information only	6.219
OSMOLALITY:	For information only	0 mOsm/kg of water
APPEARANCE:	Clear colorless to pale yellow liquid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
CELL GROWTH:	Inhibits bacterial growth without effecting cell growth. Cells have been propagated for two consecutive passages with an average of 3.9/4.0 fold increase in four days.	pass
RECOMMENDED		
STORAGE CONDITIONS:	-10 TO -20°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Michael J. Fisher

7/23/2020

Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.