

PRODUCT: Low-EDTA TE (1X), pH 8.0
10 mM Tris base, 0.1 mM EDTA

CATALOG#: 351-324-721 **Size:** 100ml

LOT#: 723834

CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
pH:	7.9 – 8.1	8.068
CONDUCTIVITY:	For information only	606 μ S
ENDOTOXIN:	Test & Report	0.048 EU/ml
APPEARANCE:	Clear, colorless liquid pass	
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass
RECOMMENDED STORAGE CONDITIONS:	15 TO 30°C	
MANUFACTURED DATE:	6/2020	
EXPIRATION DATE:	6/2023	

Michael J. Fisher

7/14/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**