

PRODUCT: DPBS with Tween 20
Dulbecco's Phosphate Buffered Saline
Without Calcium and Magnesium

CATALOG #: 114-154-151 Size: 10L

LOT #: 723829

CERTIFICATE OF ANALYSIS

	<u>RANGE</u>	<u>RESULT</u>
pH:	7.2 – 7.4	7.280
OSMOLALITY:	280 - 310 mOsm/Kg water	283
ENDOTOXIN:	Test & Report	0.022 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass

**RECOMMENDED
STORAGE CONDITIONS:** 15 – 30°C

MANUFACTURED DATE: 6/2020

EXPIRATION DATE: 6/2023

Michael J. Fisher

7/17/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**