


PRODUCT: TERRIFIC BROTH

CATALOG #: 340-071-131 Size: 1000ml

LOT #: 723822

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
pH:	7.0 – 7.4	7.140
OSMOLALITY:	For information only	525 mOsm/kg of water
APPEARANCE:	Clear, light to medium amber liquid	pass
STERILITY:	Sterilized using saturated steam under pressure	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	6/2020	
EXPIRATION DATE:	6/2022	


Quality Assurance

6/29/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**