

**PRODUCT:** WATER, Cell Culture  
**CATALOG#:** 118-162-161      Size: 20L  
**LOT#:** 723821

---

### **CERTIFICATE OF ANALYSIS**

---

This product is de-ionized by reverse osmosis and sterile filtered.

	<b><u>Specifications</u></b>	<b><u>Result</u></b>
<b>RESISTIVITY:</b>	15 – 20 Meg Ohms/cm	17.96 Meg Ohms/cm
<b>ENDOTOXIN:</b>	Less than 0.03 Eu/ml	pass
<b>MYCOPLASMA:</b>	Negative	pass
<b>APPEARANCE:</b>	Clear colorless fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
<b>CELL GROWTH:</b>	Grows mostly monolayer cells. Cells have been propagated for two consecutive passages with an average 3.0/4.0 fold increase in four days.	pass

**RECOMMENDED**

**STORAGE CONDITIONS:** 15 – 30°C  
**MANUFACTURED DATE:** 6/2020  
**EXPIRATION DATE:** 6/2022

*Michael J. Fisher*

7/09/2020

Quality Assurance

Date

FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.