


PRODUCT: Ampicillin Ready-Made Solution
100 mg/ml Ampicillin in Water

CATALOG #: 351-344-731 **Size:** 10x1ml

LOT #: 723816

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
APPEARANCE:	clear liquid	pass
PERFORMANCE:	Passed E. coli growth performance test	pass
STERILITY:	sterile filtered through 0.2 µm filter	
RECOMMENDED STORAGE CONDITIONS:	– 20°C AVOID FREEZE THAW CYCLE	
MANUFACTURED DATE:	6/2020	
EXPIRATION DATE:	6/2022	


Quality Assurance

6/18/2020
Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**