

PRODUCT: WATER, PCR QUALIFIED
DEPC Treated

CATALOG#: 351-161-741 Size: 10x1.75ml

LOT#: 723815

CERTIFICATE OF ANALYSIS

This product is ultra pure, de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Results</u>
RESISTIVITY:	> 15.11 Meg Ohms/cm	pass
APPEARANCE:	Clear, colorless liquid	pass
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass
ENDOTOXIN:	<0.06 Eu/ml	pass
MYCOPLASMA:	Negative	pass
PCR:	Not detectable for human, mouse and rat DNA contamination.	Pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	6/2020	
EXPIRATION DATE:	6/2023	

Michael J. Fisher

7/09/2020

Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.