


PRODUCT: HBSS
Hank's balanced Salt Solution
without Calcium & Magnesium

CATALOG#: 114-052-101 **Size:** 500ml

LOT#: 723780

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	6.9 – 7.3	7.076
OSMOLALITY:	266 - 294 mOsm/kg of water	278
ENDOTOXIN:	Test & Report	<0.030 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	5/2020	
EXPIRATION DATE:	5/2022	


Quality Assurance

6/9/2020

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.