

PRODUCT: Penicillin/Streptomycin Mixture
With 10,000 units/ml Penicillin
With 10,000 µg/ml Streptomycin

CATALOG#: 155-333-061 **Size:** 100ml

LOT#: 723777

CERTIFICATE OF ANALYSIS


	<u>Specifications</u>	<u>Result</u>
pH:	for information only	4.088
OSMOLALITY :	for information only	321 mOsm/Kg H ₂ O
ENDOTOXIN:	≤ 1.0 EU/ml	< 0.040 Eu/ml
MYCOPLASMA:	PCR Mycoplasma Test	pass
APPEARANCE:	Clear colorless liquid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass

RECOMMENDED

STORAGE CONDITIONS: -10 TO -20°C

MANUFACTURED DATE: 5/2020

EXPIRATION DATE: 5/2022



Quality Assurance

6/3/2020

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.