

PRODUCT: HBSS (10X)
Hank's Balanced Salt Solution
without Sodium Bicarbonate

CATALOG#: 119-065-101 Size: 500ml

LOT#: 723776

CERTIFICATE OF ANALYSIS


	<u>Range</u>	<u>Result</u>
pH:	5.5 – 6.5	5.929
ENDOTOXIN:	Test & Report	<0.030 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass

RECOMMENDED

STORAGE CONDITIONS: 15 – 30°C

MANUFACTURED DATE: 5/2020

EXPIRATION DATE: 5/2022


Quality Assurance

6/8/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**