

PRODUCT: FISCHER'S MEDIUM with L-glutamine

CATALOG #: 112-032-191 **Size:** 900ml

LOT #: 723759

CERTIFICATE OF ANALYSIS



	<u>Range</u>	<u>Result</u>
pH:	6.8 – 7.2	7.080
OSMOLALITY:	285 - 315 mOsm/kg of water	315
ENDOTOXIN:	Test & Report	<0.030 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink, to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
CELL GROWTH:	Grows mostly suspension cells - Cells have been propagated for two consecutive passages with an average 3.9/4.0 fold increase in four days.	pass

RECOMMENDED

STORAGE CONDITIONS: 2 – 8°C

MANUFACTURED DATE: 5/2020

EXPIRATION DATE: 5/2022

		6/3/2020
Quality Assurance		Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**