


PRODUCT: ACK LYSING BUFFER
CATALOG#: 118-156-101 Size: 500ml
LOT#: 723751

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.2 – 7.6	7.517
OSMOLALITY :	276 - 305 mOsm/Kg H ₂ O	302 mOsm/Kg H ₂ O
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass

**RECOMMENDED
STORAGE CONDITIONS:** 15 – 30°C
MANUFACTURED DATE: 4/2020
EXPIRATION DATE: 4/2022


Quality Assurance

5/7/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**