

PRODUCT: PBS pH 7.2
Phosphate Buffered Saline

CATALOG#: 114-056-101 Size: 500ml

LOT#: 723746

CERTIFICATE OF ANALYSIS


	<u>Range</u>	<u>Result</u>
pH:	7.1 – 7.3	7.170
OSMOLALITY:	285 - 315	302
ENDOTOXIN:	Test & Report	<0.011 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass

RECOMMENDED

STORAGE CONDITIONS: 15 – 30°C

MANUFACTURED DATE: 4/2020

EXPIRATION DATE: 4/2022



Quality Assurance

4/28/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**