


PRODUCT: DMEM Dulbecco's Modified Eagle's Medium
with 1.5 g/L Sodium Bicarbonate & with L-Glutamine

CATALOG#: 112-319-101 **Size:** 500ml

LOT#: 723744

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.0 – 7.4	7.340
OSMOLALITY:	270 - 366 mOsm/kg of water	305
ENDOTOXIN:	Test & Report	0.025 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
CELL GROWTH:	Grows mostly monolayer cells. Cells have been propagated for two consecutive passage with an average 3.0/4.0 fold increase in four days.	pass
RECOMMENDED STORAGE CONDITIONS:	2 – 8°C	
MANUFACTURED DATE:	4/2020	
EXPIRATION DATE:	4/2021	


Quality Assurance

4/27/2020

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.