


PRODUCT: PBS (1X), pH 7.4
Phosphate Buffered Saline

CATALOG#: 114-058-161 Size: 20L

LOT#: 723739

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.3 – 7.5	7.354
OSMOLALITY:	285 - 315 mOsm/Kg H ₂ O	295 mOsm/Kg H ₂ O
ENDOTOXIN:	Test & Report	<0.020 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	4/2020	
EXPIRATION DATE:	4/2022	


Quality Assurance

4/27/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**