


PRODUCT: DPBS 1X without Ca & Mg, pH 7.3-7.5
Dulbecco's Phosphate Buffered Saline

CATALOG#: 114-309-151 **Size:** 10L

LOT#: 723652

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
pH:	7.4±0.1	7.402
OSMOLALITY:	280 ± 5% mOsm/kg of water	280
APPEARANCE:	Clear colorless fluid	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	2/2020	
EXPIRATION DATE:	2/2022	



Quality Assurance

2/27/2020

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**