

**PRODUCT:** TRYPSIN (0.25%)-EDTA (0.02%) in HBSS  
without Calcium and Magnesium

**CATALOG#:** 118-093-721                      Size: 100ml

**LOT#:** 723593

---

**CERTIFICATE OF ANALYSIS**


---

	<u>Specifications</u>	<u>Result</u>
<b>pH:</b>	for information only	7.312
<b>OSMOLALITY:</b>	for information only	291 mOsm/Kg
<b>APPEARANCE:</b>	Red, to purplish-red colored fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass
<b>CELL GROWTH:</b>	Tested on monolayer cultures for cell removal and found to remove the cells in 5 minutes.	pass

**RECOMMENDED  
STORAGE CONDITIONS:** - 10 to -20°C

**MANUFACTURED DATE:** 12/2019

**EXPIRATION DATE:** 12/2021

  
Quality Assurance

1/8/2020

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**