

## Formulations

### PHOSPHATE BUFFERED SALINE

COMPONENTS	Cat.# 114-056 mg / L	Cat.# 114-058 Cat.# 130-058 mg / L	10X Cat.# 119-069 Cat.# 130-069 mg / L	Cat.# 114-154 Cat.# 130-154 mg / L
<b>Inorganic Salts</b>				
Na <sub>2</sub> HPO <sub>4</sub> (anhydrous)	726.00	795.00	7950.00	2160.00
KH <sub>2</sub> PO <sub>4</sub>	210.00	144.00	1440.00	200.00
KCl				200.00
NaCl	9000.00	9000.00	90000.00	8000.00
Tween 20®				0.5 ml/L
<b>Reference:</b> Personal Communication - Mr. George Gardiner, NIH Media Unit Tween 20® is a registered trademark of ICI Americas, Inc.				