

PRODUCT: QBSF-60 Serum-Free Medium
Stem Cell Medium
With L-Glutamine

CATALOG#: 160-204-101 Size: 500ml

LOT#: 723525

CERTIFICATE OF ANALYSIS


	<u>Range</u>	<u>Result</u>
pH:	6.8 – 7.2	7.008
OSMOLALITY:	276 - 294 mOsm/kg of water	286
ENDOTOXIN:	Test & Report	0.095 EU/ml
MYCOPLASMA:	Negative	pass
APPEARANCE:	Pink to a reddish-orange colored fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
CELL GROWTH:	Support human CD34 ⁺ cells when supplemented with the appropriate cytokines.	pass

RECOMMENDED

STORAGE CONDITIONS: 2 – 8°C

MANUFACTURED DATE: 11/2019

EXPIRATION DATE: 11/2021

 12/2/2019

Quality Assurance Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE