


PRODUCT: SeptiScan Cleanse
CATALOG#: 155-162-101 **Size:** 500ml
LOT#: 723142

CERTIFICATE OF ANALYSIS

This product is de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Result</u>
RESISTIVITY:	15 – 20 Meg Ohms/cm	17.88 Meg Ohms/cm
ENDOTOXIN:	Less than 0.03 EU/ml	pass
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
EXPIRATION DATE:	2020-05-11	


Quality Assurance

2/28/2019

Date

FOR INVESTIGATIONAL USE ONLY