



**PRODUCT:** MODIFIED FISCHER'S MEDIUM  
With Pluronic F68 and Sodium Pyruvate

**CATALOG #:** 112-121-191                      Size: 900ml

**LOT #:** 723207


**CERTIFICATE OF ANALYSIS**

	<u>Specifications</u>	<u>Result</u>
<b>PH:</b>	for information only	7.008
<b>OSMOLALITY:</b>	for information only mOsm/kg of water	317
<b>ENDOTOXIN:</b>	Test & Report	0.225 EU/ml
<b>MYCOPLASMA:</b>	Negative	pass
<b>APPEARANCE:</b>	Pink, to a reddish-orange colored fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy broth.	pass

**RECOMMENDED  
STORAGE CONDITIONS:** 2 – 8°C

**MANUFACTURED DATE:** 4/2019

**EXPIRATION DATE:** 4/2020

                      5/13/2019

---

Quality Assurance                      Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**