

**PRODUCT:** RNA GEL BUFFER 10X  
MOPS BUFFER

**CATALOG#:** 351-059-101      **Size:** 500ml

**LOT#:** 723402

---

### CERTIFICATE OF ANALYSIS


---

|                      | <u>Specification</u>    | <u>Results</u> |
|----------------------|-------------------------|----------------|
| <b>PH:</b>           | 6.9 – 7.1               | 7.019          |
| <b>CONDUCTIVITY:</b> | For information only    | 12,06 µS       |
| <b>APPEARANCE:</b>   | Clear, colorless liquid | pass           |
| <b>RNASE:</b>        | Not detectable          | pass           |

**RECOMMENDED  
STORAGE CONDITIONS:** 15 – 30°C

**MANUFACTURED DATE:** 9/2019

**EXPIRATION DATE:** 9/2022

  
Quality Assurance

10/2/2019

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**