


PRODUCT: M9 CA MEDIUM pH 7.1
CATALOG #: 340-017-101BTF Size: 500ml
LOT #: 723678

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
PH:	7.0-7.2	7.106
OSMOLALITY:	For information only	278 mOsm/kg
GROWTH PERFORMANCE ON DH11S CELLS	1.0 – 1.8	1.199
APPEARANCE:	Clear liquid	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED		
STORAGE CONDITIONS:	2 – 8°C	
MANUFACTURED DATE:	2/2020	
EXPIRATION DATE:	2/2022	


6/4/2020

Quality Assurance
Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**