



MAKING SOLUTIONS WORK

PRODUCT: Tris Glycine Transfer Buffer 10X
Transblotting Buffer

CATALOG#: 351-087-131 Size: 1000ml

LOT#: 723814

CERTIFICATE OF ANALYSIS

	<u>Specification</u>	<u>Results</u>
PH:	For information only	8.586
CONDUCTIVITY:	For information only	4.86 mS
APPEARANCE:	Clear, colorless liquid	pass
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable at 1X concentration	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Michael J. Fisher

7/22/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.**