

PRODUCT: WATER, Endotoxin Free (<0.005 EU/ml)
CATALOG#: 118-325-131 Size: 1000ml
LOT#: 723886

CERTIFICATE OF ANALYSIS

This product is de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Result</u>
RESISTIVITY:	15 – 20 Meg Ohms/cm	17.95 Meg Ohms/cm
ENDOTOXIN:	Less than 0.005 Eu/ml	pass
MYCOPLASMA:	Negative	pass
APPEARANCE:	Clear colorless fluid.	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy Broth.	pass
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass
RECOMMENDED		
STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	7/2020	
EXPIRATION DATE:	7/2022	

Chanika Johnson

 Quality Assurance

8/19/2020

_____ Date

**FOR RESEARCH USE ONLY
 NOT FOR IN-VITRO DIAGNOSTIC USE.**