



MAKING SOLUTIONS WORK

**PRODUCT:** DPBS 1X without Ca & Mg, pH 7.3-7.5  
Dulbecco's Phosphate Buffered Saline

**CATALOG#:** 114-309-151                      Size: 10L

**LOT#:** 723990

---

### CERTIFICATE OF ANALYSIS

---

	<u>Range</u>	<u>Result</u>
<b>PH:</b>	7.4±0.1	7.390
<b>OSMOLALITY:</b>	280 ± 5% mOsm/kg of water	279
<b>APPEARANCE:</b>	Clear colorless fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on Fluid Thioglycollate and Tryptic Soy broth.	pass
<b>RECOMMENDED STORAGE CONDITIONS:</b>	15 – 30°C	
<b>MANUFACTURED DATE:</b>	10/2020	
<b>EXPIRATION DATE:</b>	10/2022	

*Ta'Shemma Jackson*

10/30/2020

Quality Assurance

Date

**FOR RESEARCH USE ONLY  
NOT FOR IN-VITRO DIAGNOSTIC USE.**