



MAKING SOLUTIONS WORK

**PRODUCT:** Cryopreservation Medium  
with Fetal Bovine Serum

**CATALOG#:** 118-298-051      **Size:** 50ml

**LOT#:** 723992

---

### CERTIFICATE OF ANALYSIS

---

	<u>Range</u>	<u>Result</u>
<b>PH:</b>	For information only	7.220
<b>APPEARANCE:</b>	A reddish-orange colored fluid.	pass
<b>STERILITY:</b>	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
<b>CELL GROWTH:</b>	After the cells had been frozen and stored in liquid nitrogen for 3 days, after cells were thawed, and the viability is greater than 90%	pass

**RECOMMENDED**

**STORAGE CONDITIONS:** -10°C TO -20°C

**MANUFACTURED DATE:** 10/2020

**EXPIRATION DATE:** 10/2022

*Skin Absorption. Components of the product may be absorbed into the body through the skin. Contain DMSO. Read SDS.*

*Mycoplasma and Virus testing are performed on the appropriate raw component(s) of this medium prior to use in manufacturing.*

*Ta'Shemma Jackson*  
Quality Assurance

11/04/2020

Date

**FOR RESEARCH USE ONLY**  
**NOT FOR IN-VITRO DIAGNOSTIC USE.**